


INTUITION 2



STANDARD FEATURES

- 100% Digital
- Concealed programming socket
- Multi-memory(1-4) push button with Tone Indicator
- Direct Audio Input(DAI) capability
- Numbered Volume Control

OPTIONS

- On/Off Switch built into Battery Door
- Programmable telephone coil
- Push On replaceable earhook
- Housings available in Beige, Gray and Brown
- Battery Door Lock
- Volume Control Cover



Performance Data:

			Coupler 2cc IEC 118-7/94	Coupler MZ (711) IEC 118-0/94	Coupler 2cc ANSI S3.22-1996	Limits
SATURATION (OSPL 90)	Peak	dB SPL	134	140	134	+/- 3
	F Reference	dB SPL	127	135		+/- 3
	HF Average	dB SPL			126	+/- 3
Full-on Gain (Input: 50dB SPL)	Peak	dB	66	75	66	+/- 4
	F Reference	dB	60	70		+/- 4
	HF Average	dB			59	+/- 4
Nominal Reference Test Gain (RTG)		dB	60	60	47	
Frequency Range		Hz	200-5400			
Volume Control Range		dB	>45			
Total Harmonic Distortion at RTG:						
70 dB SPL in	500 Hz	%			4	<7
	800 Hz	%	7	7	4	<7
65 dB SPL in	1600 Hz	%			1	<3
Equivalent Input Noise Level		dB	30			<33
		dB			30	<33
Maximum Telecoil Sensitivity						
FOG; Input 10mA/m @ RTF		dB	100			
RTG; Input 31.6mA/m @ RTF		dB				
FOG; Input 31.6mA/m @ RTF		dB	110			
SPLITS @ RTF						
		dB				
HF Average		dB	103			
STS		dB	105			
		dB	1.0			
Supply Current						
	input dB SPL		60		65	
at RTG	mA		1.00		1.00	
Battery Life						
Type13 Zinc-Air(220mAh)	hrs		215		215	
Type13 Zinc-Air HighPower(200mAh)	hrs		195		195	
AGC @ 2KHz						
Attack	mS		10		10	
Release	mS		210		210	
					+/-50%	
					+/-50%	
Reference Test Frequency(RTF)						
	Hz		1600	1600	1000	

Description

The Intuition 2 BTE is a programmable 2 channel Digital Wide Dynamic Range Compression system

Highly configurable digital signal processor provides excellent versatility, with independent channel compressor characteristics

Ten-Band Gain Adjustment for precise target matching

Fully adjustable Low Frequency Slope and Corner control

Adjustable low level expansion for quieter performance

Unique Dynamic Contrast Detector "Reaches Down" to amplify quiet speech more

Look-Ahead Detection monitors signal path for loud transients and reduces gain in advance of oncoming transients, nearly eliminates signal clipping due to loud transients

Adj. Threshold Levels in each channel from 40 to 70dB


AGC-o compression limiting

Battery Type: Zinc-Air Size 13




Adjustable Low Battery Indicator and Memory Tone


Programmable with HiPro or the Micro-connect card

IEC118-0



IEC 118-7/ANSI S3.22-1996






Ten (10) Band Gain Adjustment

- Precise Target Matching
- Low Frequency Shaping
- High frequency Shaping
- Resonance Smoothing
- Feedback Notches
- Peak Shifting



Multiple Colors Available




Advanced Options

The Audina ezFIT software allows high versatility for accurate patient fittings.

- Adjustable Compression Ratio's offer a range from 1:1 to 4:1
- Adjustable and fully independant Compression and Expansion Thresholds for high flexibility
- Select from a wide range of Time Constants, MPO and Feedback Notch settings



Concealed Programming Socket




Basic Options

The Audina ezFIT software allows further fine-tuning for enhanced patient comfort.

- Adjustable 12dB/Oct. filter for High Frequency Corner control
- Adjustable Low Frequency Corner control
- Adjustable Low Slope settings of: 0dB(Flat), 6dB, 12dB and 18dB/Oct.



Multi-Memory(P1-P4) Push Button



Dynamic Contrast Detector

When two detectors are not enough, the Dynamic Contrast Detector recognizes the presence of the soft speaker and "Reaches Down" to amplify the quiet speech more.

- BASIC mode allows **very** long time constants for speech in noise
- FAST mode detects loud transients and releases quickly
- REACH mode recognizes alternating loud and soft sounds



DAI Input capability

Look- Ahead Detection

Monitors signal path for loud transients and reduces gain *in advance* of the oncoming transient- nearly eliminates signal clipping due to loud transients.



Battery Door Lock available